



This document describes the main features of the Ring BarCodeReader and computer requirements for using it.

# Ring BarCodeReader

## Main Features

1. It is intended for reading 5-character and 7-character ring barcodes.
2. The read barcode can be sent to a keyboard input buffer, which enables using the reader in another customer application without modifying it.
3. The adjustable size of the detecting area enables adjusting the readable area to the location of the barcode as much as possible.
4. It enables saving hard-to-read dosimeter images with possibilities to send them to developers for recognition algorithm improvement.



## Computer Requirements

- At least USB 2.0
- At least 4GB memory space
- Windows 7 or Windows 10 OS
- 500 MB of free disc space



25.03.2019

This specification sheet is for informational purposes only and is subject to change without notice. Windows® is trademark of Microsoft® Corporation, PDF is trademark of Adobe®, EPD®, TruDose® is trademark of Thermo Fisher Scientific.

[www.dosibase.com](http://www.dosibase.com)

Adress:

Kr. Valdemara Str. 123-40  
Rīga, LV-1013, LATVIA

Tel. +371 29214634  
Fax + 371 67283626

